



Mersa Medical
INSTRUMENTS
Manufacturer & Exporter of Quality Instruments

Manufacturer & Exporter of

Orthopedic Instruments



CONTENTS

S.NO.	DESCRIPTION	P. NOS.
1	Drill Bits & Bone Taps	1 - 2
2	Basic Instruments	3 - 7
3	Instruments For Mini Fragment	8 - 9
4	Bone Drill & Wire Instruments	10 - 13
5	Bone Holding Instruments	14 - 16
6	Bone Cutting Instruments	17 - 19
7	Bone Hammers	20
8	Hooks, Levers & Retractors	21 - 25
9	Periosteal Elevator	26
10	Instruments For Kuntscher Nails	27 - 28
11	Instruments For Hip Prosthesis	29 - 30
12	Instruments For Plaster Cutting	31
13	Instruments For Anchor Screw / Instruments For Skull Surgery & Cannulated Rigid Reamer	32
14	Distractors	33
15	Sterilization Containers & Tubes	34 - 37
16	Instruments For Interlocking Nail Removal	38



Salient Feature:

1. Pre Hardened Martensitic Stainless Steel from Japan.
2. Perfect Straightness - No Bending / Distortion Due to Heat Treatment.
3. Maximum Anti corrosive Properties.
4. Titanium Coated for Longer Drill Life.

Cat No.	Drill Bit - S.S. - Titanium Coated, Quick Coupling End & Plain Shank
TC.260.200	Size : 2.0mm Dia. X 100mm Long, Quick Coupling End
TC.260.210	Size : 2.0mm Dia. X 140mm Long, Quick Coupling End
TC.260.250	Size : 2.5mm Dia. X 115mm Long, Quick Coupling End
TC.260.270	Size : 2.7mm Dia. X 125mm Long, Quick Coupling End
TC.260.280	Size : 2.8mm Dia. X 165mm Long, Calibrated, Quick Coupling End
TC.260.320	Size : 3.2mm Dia. X 130mm Long, Quick Coupling End
TC.260.325	Size : 3.2mm Dia. X 200mm Long, Quick Coupling End
TC.260.350	Size : 3.5mm X 130mm Long, Quick Coupling End
TC.260.430	Size : 4.3mm Dia. X 221mm Long, Calibrated, Quick Coupling End

Cat No.	Cannulated Drill Bits - Titanium Coated	
TC.802.008	Size : 2.5mm Dia. X 75mm Long, Plain Shank	(For Guide Wire 0.9 MM)
TC.704.026	Size : 2.7mm Dia. X 130mm Long, Quick Coupling End	(For Guide Wire 1.25 MM)
TC.704.027	Size : 3.5mm Dia. X 130mm Long, Quick Coupling End	(For Guide Wire 1.25 MM)
TC.703.010	Size : 4.5mm Dia. X 230mm Long, Plain Shank	(For Guide Wire 2.0 MM)

Drill Bit - S.S. Quick Coupling End




Cat No.	Dia	Length
260.110	1.1mm	60mm
260.115	1.5mm	60mm
260.117	1.5mm	85mm
260.118	1.5mm	95mm
260.200	2.0mm	100mm
260.205	2.0mm	125mm
260.250	2.5mm	115mm
260.254	2.5mm	180mm
260.255	2.5mm	200mm
260.270	2.7mm	125mm
260.275	2.7mm	200mm
260.280	2.8mm	Callbrated 165mm
260.320	3.2mm	130mm
260.325	3.2mm	200mm
260.326	3.2mm	145mm
260.349	3.5mm	110mm
260.350	3.5mm	130mm
260.375	3.5mm	200mm
260.430	4.3mm	Callbrated 221mm
260.450	4.5mm	130mm
260.451	4.5mm	145mm
260.475	4.5mm	200mm
260.502	Dia. 4.5mm x 4.0mm x 310mm	

Bone Tap With T - Handle



Cat No.	Dia	Thread Length	Total Length
263.010	1.5mm	24mm	75mm
263.015	2.0mm	24mm	75mm
263.020	2.7mm	33mm	125mm
263.021	2.7mm	33mm	200mm
263.025	3.5mm	50mm	125mm
263.035	3.5mm	50mm	200mm
263.037	4.0mm	50mm	125mm
Cancellous			
263.040	4.5mm	70mm	125mm
263.045	4.5mm	70mm	175mm
263.050	6.5mm	32mm	200mm
Cancellous			

Drill Bit - Plain Shank / Jacob Chuck End



S.S. Cat No.	Dia	Length	H.S.S. Cat No.
261.611	1.1mm	60mm	-
261.615	1.5mm	60mm	-
261.620	2.0mm	100mm	-
261.622	2.0mm	200mm	-
261.625	2.5mm	115mm	262.005
261.626	2.5mm	200mm	-
261.627	2.7mm	130mm	262.010
261.629	2.7mm	200mm	-
261.632	3.2mm	130mm	262.015
261.633	3.2mm	200mm	262.020
261.635	3.5mm	130mm	262.025
261.636	3.5mm	200mm	262.030
261.645	4.5mm	130mm	262.035
261.646	4.5mm	200mm	262.040

Bone Tap Quick Coupling END



Cat No.	Dia	Thread Length	Total Length
264.010	1.5mm	24mm	50mm
264.015	2.0mm	24mm	75mm
264.020	2.7mm	33mm	100mm
264.025	3.5mm	50mm	110mm
264.035	3.5mm	50mm	175mm
264.037	4.0mm	85mm	110mm
Cancellous			
264.040	4.5mm	70mm	125mm
264.045	4.5mm	70mm	175mm
264.050	6.5mm	32mm	195mm
cancellous			
264.055	5.5mm	100mm	195mm

301.020
Hexagonal Screw Driver with Sleeve 2.5mmTip



Hexagonal Screw Driver with Self Holding Sleeve



Cat. No.	Size
301.025	For 2.5mm Tip
301.030	For 3.5mm Tip

Small Hexagonal Screw Driver



Cat. No.	Size
301.040	For 2.5mm Tip
301.050	For 3.5mm Tip

301.060
Screw Driver Shaft, Q. C. End - 2.5mmTip



301.070
Screw Driver Shaft, Q. C. End - 3.5mmTip



302.010
Small Plate Bender (Pair)



302.020
Large Plate Bender (Pair)



302.030
Plate Bending Plier (Roller Type)



Plate Bender for Reconstruction Plates



Cat. No.	Size
302.040	For 2.7mm Plates
302.050	For 3.5mm Plates
302.060	For 4.5mm Plates

**302.070
Plate Bending Press**



Depth Gauge



Cat. No.	Size
303.010	Small
303.020	Large

**304.010
Small Neutral & Loaded Drill Guide**



**304.020
Large Neutral & Loaded Drill Guide**



**304.030
DCP Drill Guide for 2.7mm Screws (2.0mm Drill Bit)**



**305.010
Counter Sink with T-Handle - 5.0mm Head
(for 1.5 & 2.0mm Screws)**



**305.020
Counter Sink with T-Handle - 6.0mm Head
(for 2.7, 3.5 & 4.0mm Screws)**



305.030

Counter Sink with T-Handle - 8.0mm Head
(for 4.5, & 6.5mm Screws)



305.038

Counter Sink - Q. C. End - 5.0mm (for 1.5 & 2.0mm Screws)



305.050

Counter Sink - Q. C. End - 6.0mm Head
(for 2.7, 3.5 & 4.0mm Screws)



305.060

Counter Sink - Q. C. End - 8.0mm Head
(for 4.5 & 6.5mm screws)



306.010

Quick Coupling Handle T - Type



306.020

Quick Coupling Handle (Fibre Handle)



306.030

Quick Coupling Adaptor for Drill Bits



306.040

Quick Coupling Tap Handle - Long , T - Type



307.040 _____
Tap Sleeve - 3.5mm



307.050 _____
Tap Sleeve - 4.5mm



307.054 _____
Double Drill Sleeve - 1.5/1.1mm



307.056 _____
Double Drill Sleeve - 2.0/1.5mm



307.058 _____
Double Drill Sleeve - 2.7/2.0mm



307.060 _____
Drill & Tap Sleeve Combined - 2.5/3.5mm



307.070 _____
Drill & Tap Sleeve Combined-3.2/4.5mm



307.080 _____
Drill & Tap Sleeve Combined-3.2/6.5mm



307.090
Insert Drill Sleeve 2.5mm



307.095
Insert Drill Sleeve 3.2mm



308.000
Pointed Drill Guide with 4.5mm Tap Sleeve & 3.2mm Insert Drill Sleeve




309.000
Screw Holding Forceps



309.002
Screw Removal Forceps



309.003
Screw Removal Forceps with Ratchet Lock



Hollow Mill



302.061
Bending Plier for Reconstruction Plates, 3.5mm & 4.5mm



Cat. No.	Size
310.010	Small - For Removal of Drill Bits
310.020	Medium - For Removal of 3.5mm Screws
310.030	Large - For Removal of 4.5 & 6.5mm Screws

Mini _____
Hexagonal Screw Driver with Sleeve (1.5mm)



Cat. No.	Length
311.001	125 mm
311.002	165 mm

311.005 _____
Mini Depth Gauge



311.006 _____
Mini Self Retaining Retractor for Fingers



311.007 _____
Mini Self Retaining Retractor for Metacarpals



311.008 _____
Mini Self Retaining Retractor for Metacarpals with Hinge



311.009 _____
Mini Neutral & Loaded Drill Guide - for 1.5mm Drill Bit



311.010 _____
Mini Plate Bender (Pair)



311.013 _____
Mini Plate Holding Forceps



311.014
Mini Plate & Bone Holding Forceps



311.015
Mini Reduction Forceps Serrated



311.016
Mini Reduction Forceps Pointed



311.017
Mini Periosteal Elevator



Mini Chisels



Cat. No.	Size
311.018	4 mm
311.019	6 mm
311.020	8 mm
311.021	10 mm
311.022	12 mm

Mini Osteotomes



Cat. No.	Size
311.023	4 mm
311.024	6 mm
311.025	8 mm
311.026	10 mm
311.027	12 mm

Mini Gouge



Cat. No.	Size
311.028	4 mm
311.029	6 mm
311.030	8 mm
311.031	10 mm
311.032	12 mm

Mini Hohmann Retractor



Cat. No.	Size
311.034	6 mm
311.035	8 mm
311.036	10 mm

320.010
Universal Bone Drill with S. S. Gears



321.020
Wire Tightner with Two Pegs



321.030
Wire Tightner Cum Twister




322.020
Guide Wire Introducer with Four Pegs



331.005
Circlage Passer



324.047
Circular Wire Cutter (Capacity up to 3.0mm)



322.030
K. Wire Bender



323.010
Wire Passer - Small



323.020
Wire Passer - Medium



323.030
Wire Passer - Large



324.010
Plier Cum Wire Bender



324.020
Plier Cum Wire Bender Cum Wire Cutter - 6"



324.030
Plier Cum Wire Bender Cum Wire Cutter - 9"



324.040
Plier Cum Wire Bender Cum Wire Cutter - 11"



324.045
Wire Cutter - Double Action, Length - 250mm



324.046
Large Pin Cutter Capacity 3.0 mm to 6.0 mm



324.050
Flat Nose Plier



325.010
Steinmann Pin Introducer with S. S. Chuck & Key



325.020
Spare S. S. Chuck & Key



325.021
Spare key for Drill Machine



325.022
Spare Chuck for Drill Machine



325.023
Spare Key for T - Handle



Steinman Pin Holder



Cat. No.	Size	Size
325.030	Small	8"
325.040	Medium	10"
325.050	Large	12"

Bone Awl with Eye




Cat. No.	Dia
326.010	2.0mm
326.011	2.5mm
326.012	3.0mm

Tripple Drill Guide



Cat. No.	Dia
328.010	For 2.0mm Wires
328.020	For 2.5mm Wires

Steinmann Pin Holder



Cat. No.	Size	Size
325.060	Extra Large	14"
325.070	Extra Large	16"
325.080	Extra Large	18"
325.090	Extra Large	20"

326.015

Reamer - Perthes, 18.5cm



Self Centering Bone Holding Forceps _____



Cat. No.	Length
340.010	150 mm
340.020	190 mm
340.030	240 mm
340.040	260 mm
340.050	280 mm

Reduction Forceps - Serrated - Ratchet Lock _____



Cat. No.	Length
341.010	140 mm
341.020	170 mm

Reduction Forceps - Serrated - Speed Lock _____



Cat. No.	Length
341.030	170 mm
341.040	200 mm
341.050	240 mm

Reduction Forceps - Pointed - Ratchet Lock _____



Cat. No.	Length
342.010	140 mm
342.020	170 mm
342.030	200 mm

Reduction Forceps - Pointed - Speed Lock _____



Cat. No.	Length
342.040	170 mm
342.050	200 mm
342.060	240 mm

343.010 _____ Patella Bone Holding Forceps - Three Prong



343.020 _____ Patella Bone Holding Forceps - Four Prong



344.010

Pelvic Reduction Forceps, Oblique, Small
Length - 190mm, with Pointed Ball Tips, Speed Lock



344.020

Pelvic Reduction Forceps, Oblique, Small
Length - 240mm, with Pointed Ball Tips, Speed Lock



344.030

Pelvic Reduction Forceps
Length - 400mm, with Pointed Ball Tips, Speed Lock



344.040

Pelvic Reduction Forceps, Small, For Screws
Length - 200mm, use with 3.5 and 4.5mm Cortical Screws



344.050

Pelvic Reduction Forceps, Medium, Length - 250mm,
use with 3.5 and 4.5mm Cortical Screws



344.060

Pelvic Reduction Forceps, Large, Length - 330mm,
Speed Lock, Fixation with 4.5mm Cortical Screws



344.070

Pelvic Reduction Forceps, with Angled Pointed
Ball Tips, Length - 250mm, Speed Lock



344.080

Pelvic Reduction Forceps, Long, with Pointed Ball Tips, Length
- 400mm, Speed Lock



344.090
Pelvic Reduction Forceps, Long, with Three Pointed Ball
Tips, Length - 400mm, Speed Lock



Lane Bone Holding Forceps with Ratchet _____



Cat. No.	Size
349.010	200mm
349.020	250mm
349.030	300mm

Lane Bone Holding Forceps without Ratchet _____




Cat. No.	Size
349.050	200mm
349.060	250mm
349.070	300mm


350.000
Burn's Bone Holding Forceps



344.100
Periarticular Reduction Forceps, Small - Jaw Opening 95mm



344.105
Periarticular Reduction Forceps, Medium-Jaw Opening 125mm



344.110
Periarticular Reduction Forceps, Large - Jaw Opening 145mm



Bone Cutting Forceps - Straight (Double Action)



Cat. No.	Size
353.010	190mm
353.020	225mm
353.030	270mm

Bone Cutting Forceps - Curved (Double Action)



Cat. No.	Size
354.010	190mm
354.020	225mm
354.030	270mm

Bone Nibbler - Straight (Double Action)



Cat. No.	Size
355.010	190mm
355.020	225mm
355.030	270mm

Bone Nibbler - Curved (Double Action)



Cat. No.	Size
356.010	190mm
356.020	225mm
356.030	270mm

Bone Nibbler - Angular (Double Action)



Cat. No.	Size
357.010	190mm
357.020	225mm
357.030	270mm

358.000 Tuder Edward Bone Cutting Forceps



359.010
Amputation Saw



360.000
Finger Saw



361.010
Bone File - Half Round - Pointed



Bone File - Flat



Cat. No.	Size
361.020	100mm
361.030	150mm

Chisels with Fibre Handle



Cat. No.	Size
362.005	5mm
362.010	10mm
362.015	15mm
362.020	20mm
362.025	25mm
362.030	30mm

Osteotome with Fibre Handle - Straight



Cat. No.	Size	Cat. No.	Size
363.005	5mm	363.016	16mm
363.008	8mm	363.018	18mm
363.010	10mm	363.020	20mm
363.012	12mm	363.025	25mm
363.014	14mm	363.030	30mm
363.015	15mm		

Osteotome with Fibre Handle - Curved



Cat. No.	Size	Cat. No.	Size
364.005	5mm	364.016	16mm
364.008	8mm	364.018	18mm
364.010	10mm	364.020	20mm
364.012	12mm	364.025	25mm
364.014	14mm	364.030	30mm

Gouge with Fibre Handle - Straight



Cat. No.	Size	Cat. No.	Size
365.005	5mm	365.016	16mm
365.008	8mm	365.018	18mm
365.010	10mm	365.020	20mm
365.012	12mm	365.025	25mm
365.014	14mm	365.030	30mm

Gouge with Fibre Handle - Curved _____



Cat. No.	Size	Cat. No.	Size
366.005	5mm	366.016	16mm
366.008	8mm	366.018	18mm
366.010	10mm	366.020	20mm
366.012	12mm	366.025	25mm
366.014	14mm	366.030	30mm

Bone Curette - Double Ended _____



Cat. No.	Size
367.010	4mm x 6mm
367.011	6mm x 9mm
367.012	9mm x 12mm

Bruns Bone Curette with Fibre Handle _____



Cat. No.	Size	Cat. No.	Size
367.023	3mm	367.028	8mm
367.024	4mm	367.030	10mm
367.025	5mm	367.032	12mm
367.026	6mm		

368.008
Bone Hammer - 100gms.



368.010
Bone Hammer - 200gms.



368.020
Bone Hammer - 350gms.



368.030
Bone Hammer - 500gms.



368.040
Bone Hammer - 700gms.



368.050
Bone Hammer - Nylon Faced



369.010
Bone Hook - Small



369.020
Bone Hook - Medium



369.030
Bone Hook - Large



370.010
Skin Hook - Single Pronged



370.020
Skin Hook - Two Pronged



370.030
Skin Hook - Three Pronged



371.000
Sharp Hook



372.010
Lane's Bone Lever - Serrated



372.020
Lane's Bone Lever - Plain - Small



372.030
Lane's Bone Lever - Plain - Large



373.000
Bristow's Retractor



374.010
Lagenback Retractor - Small



374.015
Lagenback Retractor - Medium



374.020
Lagenback Retractor - Large



375.010
Volkman's Retractor - Two Prong



375.020
Volkman's Retractor - Three Prong



375.030
Volkman's Retractor - Four Prong



375.040
Volkman's Retractor - Five Prong



376.000
Right Angle Retractor



377.000
Double Hook Retractor



379.010
Hohmann Retractor
6mm Wide, Short Narrow Point, 160mm Length



379.015
Hohmann Retractor
8mm Wide, Short Narrow Point, 160mm Length



379.020
Hohmann Retractor, 15mm Wide, 160mm Length



379.025
Hohmann Retractor
8mm Wide, Short Narrow Tip, 220mm Length



379.030

Hohmann Retractor

18mm Wide, Short Narrow Tip, 235mm Length



379.035

Hohmann Retractor

18mm Wide, Long Narrow Tip, 235mm Length



379.040

Hohmann Retractor

22mm Wide, with Wide Tip, 250mm Length



379.045

Hohmann Retractor

24mm Wide, with Long Wide Tip, 270mm Length



379.050

Hohmann Retractor

43mm Wide, Narrow Tip, 240mm Length



379.055

Hohmann Retractor

70mm Wide, Narrow Tip, 240mm Length



380.010

Self Retaining Retractor - Small - 200mm Long



380.020

Self Retaining Retractor - Medium - 250mm Long



380.030
Self Retaining Retractor - Large - 300mm Long



380.040
Self Retaining Retractor - Extra Large - 350mm Long



380.050
Gelpi Retractor - Small - 175mm



380.060
Gelpi Retractor - Medium - 245mm



380.070
Gelpi Retractor - Large - 315mm



Periosteal Elevator - Straight



Cat. No.	Size	Cat. No.	Size
381.006	6mm	381.014	14mm
381.008	8mm	381.016	16mm
381.010	10mm	381.020	20mm
381.012	12mm	381.025	25mm

Periosteal Elevator - Curved

Cat. No.	Size	Cat. No.	Size
382.006	6mm	382.014	14mm
382.008	8mm	382.016	16mm
382.010	10mm	382.020	20mm
382.012	12mm	382.025	25mm



381.201

Periosteal Elevator - 6.0mm Width, Straight Blade, Straight Edge



381.202

Periosteal Elevator - 13.0mm Width, Straight Blade, Straight Edge




381.203

Periosteal Elevator - 14.0mm Width, Straight Blade, Round Edge



382.201

Periosteal Elevator - 3.0mm Width, Curved Blade, Straight Edge



382.202

Periosteal Elevator - 6.0mm Width, Curved Blade, Straight Edge



382.203

Periosteal Elevator - 6.0mm Width, Curved Blade, Round Edge



382.204

Periosteal Elevator - 13.0mm Width, Curved Blade, Straight Edge




382.205

Periosteal Elevator - 14.0mm Width, Curved Blade, Round Edge



Bone Spreader _____



Cat. No.	Size
383.014	140mm
383.021	210mm
383.027	270mm

384.004 _____
K. Nail Setting Device



384.008 _____
K. Nail Measuring Gauge



Guide Wire for K. Nail _____



Cat. No.	Size
384.012	2.5mm for 6 & 7mm Nails
384.014	3.0mm for 8 & 9mm Nails
384.016	4.0mm for 10 & 11mm Nails
384.018	4.5mm for 12 & 13mm Nails

384.002 _____
K. Nail Driver




384.006 _____
K. Nail Extractor with Two Hooks



384.010 _____
Diamond Pointed Awl



K. Nail Reamer _____



Cat. No.	Size	Cat. No.	Size
384.020	6mm	384.030	11mm
384.022	7mm	384.032	12mm
384.024	8mm	384.034	13mm
384.026	9mm	384.036	14mm
384.028	10mm		


384.040
Spare Hook for K. Nail Extractor



384.042
Slide Hammer for K- Nail Extractor Set



384.050
Locking Plier - Large



384.060
Extraction Rod with Hammer for Locking Plier



386.010
Rasp for Austin Moore - Standard Stem



386.011
Rasp for Austin Moore - Narrow Stem



386.020
Rasp for Thompson - Standard Stem



386.021
Rasp for Thompson - Narrow Stem



386.030
Murphy Lane Bone Skid



386.040
Moore Hollow Chisel



386.041
Moore Hollow Chisel with Fibre Handle - Standard



386.042
Moore Hollow Chisel with Fibre Handle - Extra Large



386.050
Judet Auger Extractor



Measuring Gauge for Prosthesis



Cat. No.	Size
386.061	Set of Ten (38 to 56mm)
386.062	Set of Ten (37 to 55mm)

386.070
Impactor - Nylon Faced



386.071
Impactor - Nylon Faced, Large Tip



386.080
Extractor with Two Hooks



386.090
Spare Hook for Extractor



388.005
Plaster Saw Heavy Handle



388.010
Plaster Saw Engel



388.015
Plaster Saw Bergman



388.020
Plaster Spreader



388.025
Plaster Shear - Bohler's



388.030
Plaster Shear - Stille's



388.035
Plaster Shear - Gay's



388.040
Bandage Cutting Scissor



388.045
Plaster Bender



388.050
Plaster Opener - Daw's



802.050
Anchor Screw Driver



802.051
Anchor Guide Sleeve



Anchor Pointer



Cat No.	Size
802.052	3.0mm
802.053	5.0mm

802.054
Suture Puller With Loop



Instruments for Skull Surgery

Gigli Saw Wire
50cm long



Cat. No.	Size
394.005	Ordinary
394.010	Fine (Imported)

394.020
Gigli Saw Introducer



394.015
Gigli Saw Handle - Pair



Cannulated Rigid Reamer



Cat. No.	Size	Cat. No.	Size
708.201	6.0mm	708.209	11.0mm
708.202	7.0mm	708.210	11.5mm
708.203	8.0mm	708.211	12.0mm
708.204	8.5mm	708.212	12.5mm
708.205	9.0mm	708.213	13.0mm
708.206	9.5mm	708.214	13.5mm
708.207	10.0mm	708.215	14.0mm
708.208	10.5mm		

396.010
Small Distractor for Tibia



396.015
Large Distractor for Tibia & Femur



396.101
Distractor Single hole 50mm



396.102
Distractor Double Hole 50mm



Distractor Double Hole



Cat. No.	Size
396.201	4.0mm x 75mm
396.202	4.0mm x 100mm
396.203	4.0mm x 150mm
396.301	6.0mm x 150mm
396.302	6.0mm x 200mm
396.303	6.0mm x 250mm
396.304	6.0mm x 300mm
396.305	6.0mm x 350mm

400.005
Instrument Box with Two Trays Length 300mm



400.010
Small Fragment Implant Box with Two Trays for Instruments and One Tray for Screw with Screw Holding Forceps.



400.015
Instrument Box with Two Trays - Length 450mm



400.020
Large Screw Box With Screw Tray & Screw Holding Forceps for 4.5mm Cortical, 4.5mm Malleolar & 6.5mm Cancellous Screws



400.025
Small Screw Box With Screw Tray & Screw Holding Forceps for 3.5mm Cortical, 3.5 & 4.0mm Cancellous Screws



Screw Box for Interlocking Bolts With
Screw Tray & Screw Holding Forceps



Cat. No.	Size
400.030	for 4.9 mm Locking Bolts
400.031	for 3.9 mm Locking Bolts
400.032	for 3.4 mm Locking Bolts

400.033
Combined Screw Box for 3.4, 3.9, 4.9mm Interlocking Bolts
With Screw Tray & Screw Holding Forceps



400.035
Screw Box Combine with Large & Small Screw Tray & Screw
Holding Forceps for 3.5, 4.0, 4.5, 6.5mm Screws



400.040
DHS/DCS Implants Box with Three Tray for DHS/DCS
Screws & Plates



400.045
Sterilization Box for Flexible Reamers with Detachable Heads




400.050
Sterilization Box for Fixed Head Flexible Reamers with Two
Trays S.S. L-210mm W-155mm H-60mm



400.055
Instrument Box for PFN & Re-con Nails



400.056
Graphic Box for Delta I.L. Set (Tibia & Femur)



400.059
Screw Box for 4.0mm Cannulated Cancellous Screw - with Screw tray & Forceps.




400.060
Screw Box for 7.0mm Cannulated Cancellous Screws with Screw Tray, Instrument Tray & Screw Holding Forceps



400.061
Screw Box - for 3.5mm & 4.0mm Locking Head Screw



400.062
Screw Box - for 5.0mm Locking Head Screw



400.073
Box for PFN/TFN Bolts & I.L. Bolts with Screw Tray & Screw Holding Forceps.




Sterilizing Tubes
For Sterilization of Wires, Nails, Pins & Drill Bits



Cat. No.	Size
401.002	150mm
401.004	225mm
401.006	300mm
401.008	400mm
401.010	500mm

400.070
Sterilization Box For Interlocking Nails Removal S.S.



723.101
Extractor Rod




723.102
Handle For Extractor



723.103
Hammer For Extractor Rod



723.105
T-Wrench 10mm Cannulated



Conical Bolt




Conical Extraction Reamer (For Removal of I.L. Nails)



Cat. No.	Size
723.120	20mm Tapped Thread TPI
723.122	22mm Tapped Thread TPI
723.124	24mm Tapped Thread TPI
723.126	26mm A Tapped Thread TPI
723.128	26mm B Tapped Thread TPI
723.116	16mm Straight Thread TPI

Cat. No.	Size
310.360	9,10,11mm Nails
310.361	12,13,14mm Nails

IS.052
Sterilization Box for Interlocking Nails Removal S.S.
(Including Instruments)





Roras Road, New Adalat Grah
Sialkot-51310 Pakistan

Cell: +92-321 6115424

Email: info@mersamedical.com

Email: shahzad@mersamedical.com

Web: www.mersamedical.com